

Model	Go Chair®
Part Number	GoChair
Weight Capacity	300 lbs.
Maximum Speed ¹	Up to 3.7 mph
Ground Clearance ³	1.8" at motor
Turning Radius ³	25.5"
Overall Length ³	33" with front rigging
Overall Width ³	25.5"
Front Wheels	6" solid casters
Drive Wheels	9" solid
Rear Casters	3" solid anti-tips
Suspension	Limited
Range Per Charge ^{1,2}	Up to 13.2 miles (200 lbs.) Up to 8.7 miles (300 lbs.)
Heaviest Piece When Disassembled ³	36 lbs. (rear section)
Standard Seating	Type: Compact, foldable Weight: 30 lbs. Material: Black vinyl Dimensions: Width: 18", Depth: 17"
Standard Electronics	40 amp, PG GC3 (non-programmable)
Battery Requirements ⁴	(Two) 18 AH, 12 V, deep-cycle Weight: 14 lbs. each
Battery Charger	Off-board, 2 amp
Warranty	Frame: 5-year limited Drivetrain: 13-month limited Electronics: 13-month limited Batteries: 13-month limited

HCPCS Code: K0899

POWER MOBILITY DEVICE, NOT CODED BY DME PDAC
OR DOES NOT MEET CRITERIA

- ¹ Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tire condition. This specification can be subject to a variance of (+ or -) 10%.
- ² Tested in accordance with ANSI/RESNA, MC Vol2, Section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
- ³ Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.
- ⁴ Battery weight may vary based on battery manufacturer.

PRIDE[®]
MOBILITY
LIVE YOUR BEST[®]

(US) 800-800-8586 • (Canada) 888-570-1113

pridemobility.com |    

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.

GOCHAIR_061819

GOCHAIR[®]

**GET
OUT
THERE!**



LIVE YOUR BEST[®]

pridemobility.com |    

GET OUT THERE

WITH THE #1 PORTABLE POWER CHAIR!



Disassemble quickly for transport or storage.



Store items close at hand with convenient under seat storage.



Personalize your Go Chair® with one or more optional accessories.

The Go Chair® is #1 for portability! With convenient feather-touch disassembly, the Go Chair allows you to experience lightweight travel and independence on the go. Enjoy dual swing-away storage bins and a compact turning radius for outstanding maneuverability. Available in an array of contemporary colors, it's time to get out there with the Go Chair!

Additional Features

- 18" x 17", 5-point swivel seat with latch release
- Compact turning radius
- Dual swing-away, under seat storage bins
- Height and width adjustable armrests
- Large footplate
- One hand feather-touch disassembly
- 60" lap belt
- The heaviest piece is only 36 lbs.

Available Colors



Accessories*

- Armrest mount cane clip
- Cane/crutch holder (3" single, 3" double)
- Cup holder
- Lap belt (50", 70")
- Oxygen tank holder
- Rear mount bracket
- Saddlebag
- Safety flag
- Weather cover
- XLR USB charger
- Additional battery pack

Military Patches



* Only one rear accessory can be used at a time.